

Medical Device Industry Overview

E-Learning course PDF



Copyright

Copyright © 2017 Honeywell International Inc. All rights reserved. Honeywell and its product names are among the trademarks and/or service marks owned by Honeywell International, Inc., or its subsidiaries. All other product names mentioned herein are trademarks or registered trademarks of their respective owners, who may or may not be affiliated with, connected to, or sponsored by Honeywell. Reference to any product, process, publication, service, or offering of any third party by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply the endorsement or recommendation of such by Honeywell International Inc.

Published By

Honeywell Safety and Productively Solutions 703 Rodi Road, Pittsburgh, PA 15235 (412) 829-8145 fax (412) 829-0972 http://www.vocollectvoice.com

Confidentiality

For use by employees, partners, and customers of Honeywell—this documentation provides information for customers who are using Honeywell products.

Warning and Disclaimer

Honeywell International Inc. ("HII") reserves the right to make changes in specifications and other information contained in this document without prior notice, and the reader should in all cases consult HII to determine whether any such changes have been made.

The information in this publication does not represent a commitment on the part of HII. HII shall not be liable for technical or editorial errors or omissions contained herein; nor for incidental or consequential damages resulting from the furnishing, performance, or use of this material.

This document contains proprietary information that is protected by copyright. All rights are reserved. No part of this document may be photocopied, reproduced, or translated into another language without the prior written consent of HII.



Contents

Instructions	
Introduction	
Welcome	
Course Objective	6
Medical Devices move to the home	
Hemodialysis	
Peritoneal dialysis	
Growth in portable- and home-use medical devices	
Internet of things- connected medical equipment	10
Target/focus medical applications	11
Kidney dialysis machines	
Respiratory	
Diagnostic/analytical equipment	
Hospital hardware	15
Dental equipment	
Surgical instruments	17
Patient monitoring systems	
Infusion pumps	
Medical applications for Honeywell sensors-Activity	20
Course Summary	2 1



Instructions




Course instructions



This course has audio and video. Please put on your headset and make sure the volume is set to a comfortable level. You can hide or display the Transcript and Topics at any time during the course. A PDF version of this course is available to download or print from the Print section. For any questions or issues, please email Training. The duration of this course is 20 minutes. Select Start to proceed with the session.




Introduction



Hello and welcome to the Medical Device Industry Overview course. This course introduces trends and technological changes that influence the medical device industry and the role of Internet of Things, or IOT in the evolution of this industry.



Course Objective




After completing this course, you will be able to:

- Describe how medical device industry trends and technological advances are relevant to sensors,
- Identify how this is shifting to consumer-like behaviors,
- Describe the impact of Internet of Things, or IOT, on medical equipment, and
- List medical applications relevant to Honeywell 's sensor portfolio



Medical Devices move to the home



Select each label shown on the screen to learn more.

Hemodialysis

Treatment models that were once limited to facilities are now available in other settings, thanks to technology. As a result, medical devices are becoming more portable and accessible for everyone. For example, look at dialysis. Dialysis removes impurities from blood when a patient's kidneys are unable to perform this critical, life-sustaining function.


Would you prefer travelling to a facility multiple times every week for treatment or have treatment performed in your home?

Would you prefer needles or a surgically implanted catheter?

Today, critical, life sustaining treatments like dialysis are moving to the home environment. Rather than traveling to a hospital or clinic three to four times a week for hemodialysis, where a needle and tubing connect the patient to a machine that cleans and returns the blood to the patient, many people are now able to have dialysis performed in the comfort of their own homes. This is called peritoneal dialysis.



Peritoneal dialysis



Peritoneal dialysis doesn't require blood removal or return. The patient uses the catheter to administer fluid into the abdomen. This fluid removes impurities from the blood and then drains into a collection bag. Peritoneal dialysis eliminates the need to travel to a clinic, use needles, or obtain nurse support.



Growth in portable- and home-use medical devices



Diagnostics, monitoring, and even treatments are moving to the home environment where possible, as it reduces healthcare costs and is more convenient for the patient.

For years, devices like glucometers, heart rate, and blood pressure monitors have been available in the home. More advanced versions of these devices transmit vital sign data from a patient's home to clinicians, enabling real-time monitoring of a patient's health.


More advanced treatment devices, like dialysis machines, sleep apnea machines, infusion pumps, and oxygen concentrators are available in the home environment.

All of these devices utilize some kind of sensor. For example, blood pressure monitors use pressure sensors. Peritoneal dialysis machines use pressure sensors and temperature sensors. Infusion pumps use force sensors.

Honeywell is poised to take advantage of increased sales volume in the sensor market from the growth in home-based medical devices.



Internet of things- connected medical equipment



Medical devices are becoming more connected.

The proliferation of medical devices with embedded sensors, combined with existing connectivity infrastructure, has fostered a dramatic increase in health-related applications and real-time data monitoring. This, in turn, fuels the development of additional connected devices.


Honeywell has various sensor offerings, many with industry-leading, small-size, and power-consumption features, making them ideally suited to mobile and connected devices. Let's take a look.



Target/focus medical applications

Target/focus medical applications

Select each medical application to learn more about it.



Medical devices are used in a wide span of medical applications. Honeywell offers a variety of sensors and switches focused on eight medical applications with high win rates. Note that these aren't the only areas we can sell into, but these are the focus areas. We do sell into some devices that are outside of the eight target/focus medical applications.

Select each medical application shown on the screen to learn more about it.



Kidney dialysis machines

Target/focus medical applications

Select each medical application to learn more about it.

Respiratory

Diagnostic/analytical equipment

Hospital hardware

Dental equipment

Surgical instruments

Patient monitoring systems

Infusion pumps

- Replaces some kidney functions
- Removes bloodstream waste and fluid through
 - Diffusion
 - Osmosis of solutes and fluid
- · Examples of medical equipment:
 - Hemodialysis machines
 - · Peritoneal dialysis machines



Kidney dialysis machine

A kidney dialysis machine is used to replace some kidney functions by removing waste and fluid from the bloodstream via diffusion and osmosis of solutes and fluid across a semi-permeable dialysis membrane. Peritoneal dialysis and hemodialysis machines are examples of medical equipment in this category.

Respiratory

Target/focus medical applications

Select each medical application to learn more about it.

Respiratory

Diagnostic/analytical equipment

Hospital hardware

Dental equipment

Surgical instruments

Patient monitoring systems

Infusion pumps

- Detect/corrects respiratory interruptions
- · Examples of medical equipment:
 - · Anesthesia delivery systems
 - Ventilator
 - Spirometer
 - Continuous Positive Airway Pressure (CPAP) machine



Respiratory applications

The respiratory medical application represents a set of medical equipment responsible for uninterrupted respiration, that is, inhale of oxygen and exhale of carbon dioxide from the human body. Sensors are used to detect any interruptions in the airflow to and from the body.

Some examples of medical equipment in this category include anesthesia delivery systems, ventilators, spirometers to measure lung function, and Continuous Positive Airway Pressure, or CPAP machines for sleep apnea.



Diagnostic/analytical equipment

Target/focus medical applications

Select each medical application to learn more about it.

Respiratory

Diagnostic/analytical equipment

Hospital hardware

Dental equipment

Surgical instruments

Patient monitoring systems

Infusion pumps

- · Provide data to:
 - Diagnose diseases
 - · Facilitate decision making
 - Monitor health
- · Examples of medical equipment:
 - Blood analyzer
 - Chemistry analyzer



Diagnostic/analytical equipment

Diagnostic and analytical equipment is used to diagnose diseases, make treatment decisions, and to monitor patients.

Some examples of medical equipment in this category include blood analyzers and chemistry analyzers.

Hospital hardware

Target/focus medical applications

Select each medical application to learn more about it.

Respiratory

Diagnostic/analytical equipment

Hospital hardware

Dental equipment

Surgical instruments

Patient monitoring systems

Infusion pumps

- Facilitate smooth hospital operations
- · Examples of medical equipment:
 - Sterilizers
 - Blood storage refrigerators
 - Autoclaves
 - Incubators
 - Hospital beds
 - Medication dispensing cabinets



Hospital hardware

Hospital hardware facilitates smooth hospital operations and patient care.

Some examples of medical equipment in this category include sterilizers, blood storage refrigerators, autoclaves, incubators, hospital beds, medication dispensing cabinets etc.

Dental equipment

Target/focus medical applications

Select each medical application to learn more about it.

Respiratory

Diagnostic/analytical equipment

Hospital hardware

Dental equipment

Surgical instruments

Patient monitoring systems

Infusion pumps

- Facilitate smooth dental operations
- Regulate motion control and positioning
- · Examples of medical equipment:
 - Dental imaging systems
 - Dental chairs
 - Pressure-operated dental instrument—drills, water sprays, and air blasters



Dental equipment

Dental equipment uses sensors for position control of dental chairs and pressure monitoring of fluid or air.

Some examples of medical equipment in this category include dental imaging systems, dental chairs, and pressure-operated dental instruments including drills, water sprays, and air blasters.

Surgical instruments

Target/focus medical applications

Select each medical application to learn more about it.

Respiratory

Diagnostic/analytical equipment

Hospital hardware

Dental equipment

Surgical instruments

Patient monitoring systems

Infusion pumps

- · Ensure patient safety
 - Regulate movement of air and gas
 - Maintain precise pressure and force levels
- Examples of medical equipment:
 - o Orthopedic bone drills
 - Endoscopes



Surgical instruments

Surgical instruments use sensors to regulate the movement of air, gas, or fluids to maintain precise levels of pressure and force in surgical instruments.

Some examples of medical equipment in this category include orthopedic bone drills, endoscopes, and surgical fluid management systems.

Patient monitoring systems

Target/focus medical applications

Select each medical application to learn more about it.

Kidney dialysis machines

Respiratory

Diagnostic/analytical equipment

Hospital hardware

Dental equipment

Surgical instruments

Patient monitoring systems

Infusion pumps

- · Facilitate health monitoring
 - · Regulate airflow pressure
 - · Monitor blood pressure
 - Monitor temperature
- · Examples of medical equipment:
 - Blood pressure monitors
 - o Blood glucose monitors
 - · Respiratory monitors
 - · Temperature monitors



Patient monitoring &systems

Patient monitoring equipment is used in operating rooms, emergency rooms, intensive care units, and patients' homes to monitor and display vital signs of patients' health conditions, such as blood pressure, blood glucose, respiration, temperature etc. This data enables doctors and medical practitioners to make treatment decisions.

Some examples of medical equipment in this category include patient monitoring, blood pressure and blood glucose monitoring, respiratory monitoring, and temperature monitoring.

Infusion pumps

Target/focus medical applications

Select each medical application to learn more about it.

Respiratory

Diagnostic/analytical equipment

Hospital hardware

Dental equipment

Surgical instruments

Patient monitoring systems

Infusion pumps

- · Ensure delivery of fluids
 - Uninterrupted
 - Accurate
 - Controlled
- · Detect occlusion or blockage
- · Examples of infusion pumps:
 - High-volume infusion pumps
 - Insulin pumps
 - · Syringe pumps




Infusion pumps

Infusion pumps are used to deliver fluids, such as nutrients or medications, into a patient's body accurately and in a controlled manner. Sensors are used to monitor for occlusion or blockage preventing the delivery of fluids to the patient.

Some examples of medical equipment in this category include high-volume infusion pumps, insulin pumps, and syringe pumps.

Let's see how much have you retained from the information shared on this screen.

Medical applications for Honeywell sensors-Activity






In this activity, labels with different medical equipment names will appear on the screen. Each medical equipment is an example of one of the focus medical applications. Drag each medical equipment label on to the correct medical application name listed on the screen. After making a drop, select Proceed to view the next medical application label and its drop options. There are 16 opportunities and you need at least 12 correct matches to complete this activity. There are unlimited attempts and the duration of



each attempt is four minutes. Note that after each drop, you will need to select Proceed on the screen to review the next drop and its drop options. Select Start to begin.

Course Summary



In this course, you learned that:

Medical devices have become smaller, portable, and personalized offering people convenience of operating and using them in the home environment.

Small portable devices mean increased opportunities in the sensor sales market.

Internet of things, or IOT, has made gathering medical data for current and future analysis an effortless activity through digitization and interconnectivity of medical devices.

There are eight medical applications aligned to Honeywell product offerings including: kidney dialysis machines, respiratory, infusion pumps, diagnostic or analytical equipment, patient monitoring systems, hospital hardware, surgical equipment, and dental equipment.

